
REGIONAL ANAESTHESIA
Peripheral blocks

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Value Life

ENGLISH - pp. 1-20

- Electric cable with three poles plug (1a), Electrode lead (1b) and Needle lead (1c).
- Display
 - Current or charge;
 - Sound tones (enabled or disabled)
 - Precision feature for fine adjustments (enabled);
 - Side control (enabled);
 - Units of measurement selectable options;
 - Frequency selectable options;
 - Impulse duration selectable options;
 - Unit of measurement selected;
 - Exhausted battery indicator
- Rotating knob for adjustments of current or charge supplied
- Set up button or Modification of impulse duration button
- Side control for adjustments of current or charge supplied
- Safety button to zero immediately the current or the charge supplied
- ON/OFF button
- Battery panel (in the back lower side of the device)

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INDEX

1. OVERVIEW	2
1.1 Intended use	2
1.2 Warnings	2
2. FUNCTIONING OF THE STIMULATOR	4
2.1 Set up of the stimulator	4
2.1.1 Check list of the parts	4
2.1.2 Battery installation / change	4
2.1.3 Cables connection	4
2.2 Stimulator's switch on / off	4
2.3 Set up modification	5
2.4 Current / charge selection	5
2.5 Display	6
2.6 Side Control	7
2.7 Safety features	7
2.7.1 Current self limitation at switching on	7
2.7.2 Sound tone	7
2.7.3 Set / supplied current	8
2.7.4 Disconnection alarm	8
2.7.5 "Safety" button	8
2.7.6 "Precision" feature	8
2.7.7 Exhausted battery	8
2.7.8 Self switching off	9
3. INSTRUCTIONS FOR USE	9
3.1 Use of the device	9
3.2 Check of the functioning	10
3.3 Conditions of use, cleaning and maintenance	10
3.4 Accessories	11
4. TECHNICAL DATA	12
4.1 General information	12
4.2 Output waveform	12
4.3 Precision	13
4.4 Classification	13
4.5 Liability	13
4.6 Waste disposal and conformity to the WEEE directive	14
4.7 Transport and storage	14
4.8 Batch and year of manufacturing	14
4.9 Electromagnetic compatibility	15
5. REGULATORY ISSUES	19
5.1 Applied standards	19
5.2 Declaration of conformity	19
6. WARRANTY	20

I. OVERVIEW

This device is manufactured according to current International safety standards. Please carefully read the following warnings, instructions and safety information before using the device as recommended by the symbol  printed on the device.

1.1 Intended use

Plexygon is an electronic medical device aimed to provide the anaesthetist, practising peripheral anaesthesia or pain therapy, with an instrument allowing an easier and more precise localisation of the nerve. Nerve location is achieved by generating an appropriate electrical stimulation, through a dedicated needle-electrode, then inducing a muscular reaction on the afferences of the nerve reached by the stimulation. The PLEXYGON stimulator has been designed to obtain the best performances in term of reliability, user friendliness, efficacy and safety.

1.2 Warnings

The warning notes are marked along the following paragraphs by the symbol 

Check of the functioning

- If Plexygon is not working as described in this manual, don't use the device and call for technical assistance.
- The user must make sure that the device is safe and integral before using it. Carry out verifications as specified (see paragraph 3.2).
- If at the time of switching on the display is any different from what is described, the stimulator is not suitable for the use.
- If the symbol of battery is blinking, this means that the battery is not charged enough and must be changed.

Use of the stimulator

- The stimulator must be used by qualified personnel and only after completely reading the enclosed manual of instructions for use. These personnel must also have a sound knowledge of anatomy and neuro-physiology, and knowledge of the protocols in force for the execution of peripheral blocks (see paragraph 3.1).
- The injection of anaesthetic drugs causes a variation of the nerve response to the stimulation; therefore after the injection it will not be possible to maintain the stimulation at the same levels of current.-Always test for aspiration of blood before any injection.
- Using lower charge allows further precision but implies higher risks of direct contact with the nerve.
- If the sound tone has been deactivated, check the display to make sure that the current is effectively supplied to the patient (see par. 2.7.2).
- Since the device is not supposed to be sterile, do not bring it in a sterile field

- nor get it in touch with sterile personnel.
- The use of the stimulator can interfere with the functioning of Pacemakers or other stimulators implanted or anyhow connected to the patient. Therefore don't use the stimulator without consulting a specialist of these devices.
 - Mobile, radio frequency communication devices (including cables for external antennas and antennas) must not be used closer than 30 cm (12 inches) to any part of the device, including the cables specified by the manufacturer. Otherwise, the performance of the device may be affected.
 - Use of this device near or over (above or below) other devices is prohibited as it may cause improper functioning. If this type of use is necessary, the device in question and the others involved must be closely monitored to verify their correct functioning.
 - The simultaneous connection of the patient to a surgical high frequency device can cause burns in correspondence to the electrodes of the stimulator, or anyway damages to the stimulator itself.
 - Don't use the stimulator in presence of flammable gases.
 - Verify that the cables don't touch with other cables or directly with the patient.
 - Don't use the stimulator in case of any malfunctioning, leakage of battery fluids, accidental pouring of liquids on the stimulator.
 - Don't use non-original accessories. The use of non-original accessories can modify the emission of the electrical impulse, increase electromagnetic emissions and/or decrease the electromagnetic immunity of the device and consequently affect the safety features of the stimulator itself.
 - Do not act near microwave and/or short waves therapy devices.
 - High risk of cardiac fibrillation due to application of electrodes in the thoracic area.
 - Check the compatibility between the device and its accessories and handle with particular care.
 - The stimulation should not be applied across or through the head, directly on the eyes, covering the mouth, on the front of the neck (especially the carotid sinus), or from electrodes placed on the chest and the upper back or crossing over the heart.
 - Do not touch the battery contacts and the patient at the same time.
 - To reduce the risk of liquid leakage from the battery, remove the battery from the device when not in use for a long time.

Cleaning and maintenance

- The stimulator must be cleaned by a soft cloth dampen with natural water; it's possible to use also degreasing or sterilizing substances. The usage of ethylic alcohol is NOT allowed. Don't use sprays nor immerse the stimulator in any liquid.
- A spilling or a stagnation of liquid or humidity on the device can cause a malfunction of the stimulator, therefore it is important to avoid any spills of liquid on the stimulator. In case of accidental spilling mop carefully any eventual residue.
- For the device disposal please refer to the instructions in paragraph 4.6.

2. FUNCTIONING OF THE STIMULATOR

2.1 Set up of the stimulator

2.1.1 Check list of the parts

Verify that inside the case, further to the stimulator, the following accessories are included:

- 1 cable for connecting the stimulator to the needle.
- 1 alkaline battery 9V - 500 mAh (type 6LR61);

2.1.2 Battery installation / change

Open the panel in the back lower side of the stimulator by gently pushing on the lever of the panel itself and pulling upwards.

Then connect an alkaline battery 9V-500mAh (type 6LR61) to the clip inside the battery room, caring of the proper connection of the positive and negative poles. Insert the battery in its room and close the panel with a procedure opposite to the opening one.

2.1.3 Cables connection

Once the battery has been inserted and the back panel closed, set up the double lead cable by inserting the three poles plug in the stimulator's front socket. Thereafter connect the BLACK lead in to the needle (negative electrode; black plug) and to the ECG electrode RED lead on the patient electrode.

2.2 Stimulator's switch on / off

To switch the stimulator on, press the ON / OFF button on the front panel.

By that the display will show the default set up: the current associated to the supplied impulse (0,00 mA at the switching on), the frequency (2 Hz) and the duration of the stimulus itself (300 μ S).

Before proceeding further, it's recommended to verify the proper functioning of the stimulator according to the modalities described in par. 3.2 (Check of the functioning).

To switch the stimulator off, keep the ON/OFF button pressed for 2 sec.

2.3 Set up modification

µS
MODE

To modify the set up of the stimulator; keep the $\mu\text{S}/\text{MODE}$ button pressed for 2 sec. Once entered in the menu, the modifiable parameter is pointed; then it's possible to change the setting by rotating the knob clockwise or anticlockwise. Pressing again the $\mu\text{S}/\text{MODE}$ button, the selection is confirmed and the next parameter is selected.

The sequence of the selection is the following:

Frequency: Selectable options: 2 - 1 - 4 Hz
 Duration: Selectable options: 300 - 50 - 100 μS
 Measure unit: Selectable options: mA - nC
 Sound: Selectable options: on - off

A further push of the MS/MODE button allows to go out from the options menu and return to the normal functioning.

During the functioning of the stimulator; it is possible to modify the duration of the impulse keeping constant the charge supplied (expressed in nC) with a pressure and an immediate release of the $\mu\text{S}/\text{MODE}$ button. If the visualisation is in mA, the current will change in order to keep constant the charge supplied (for example: duration 300 μS - current 0.5 mA; modified duration to 100 μS - current to 1.5 mA; modified duration to 50 μS - current to 3.0 mA)

2.4 Current / charge selection

The knob on the upper right side of the front panel allows adjustment of the current or the charge supplied by the stimulator at any impulse.

At each rotation step, the selected current will change by 0.1 mA or 0.02 mA respectively if its value is bigger or lower than 0.5 mA; turning the knob clockwise the current increases, turning anticlockwise the current decreases. During the fine adjustments (under 0,5 mA) the symbol "P" (Precision) appears on the display.

NOTE: If the visualization in "nC" has been selected, the charge variation step is equal to the product of the current multiplied by the duration of the impulse, corresponding then to the values indicated in the table here below:

Table I: Variation steps of the charge Q (nC) for different durations of the impulse

	300 μS	100 μS	50 μS
$Q \geq 150 \text{ nC}$	30 nC	-	-
$Q < 150 \text{ nC}$	6 nC	-	-
$Q \geq 75 \text{ nC}$	-	10 nC	-
$Q < 75 \text{ nC}$	-	2 nC	-
$Q \geq 25 \text{ nC}$	-	-	5 nC
$Q < 25 \text{ nC}$	-	-	1 nC

The range of selectable current is:

from 0 to 4 mA for impulse duration of 300 μs (or 0 to 1200 nC),

from 0 to 5 mA for impulse duration of 100 μs (or 0 to 500 nC),


from 0 to 6 mA for impulse duration of 50 μs (or 0 to 300 nC).

Each delivered impulse is accompanied by a sound.

The sound has a tone increase each 0,1 mA between 1,0 and 0,5 mA.

For safety reasons it is not possible, for any current or duration, to set a current higher than 2,0 mA (or an equivalent charge) before the needle has been introduced in the patient (see par. 2.7.1)

2.5 Display

The display is used either for the visualization of the current or charge associated to each impulse, and for showing the selected parameters of functioning (duration of the impulse, frequency, measure unit of the visualization, active selector, and eventual battery alarm).

The visualized figure corresponds to the impulse actually supplied by the stimulator; apart from when the knob

is turned. In this case the figure shown corresponds to the selected figure during the adjustment and for the following four seconds.

If the supplied current differs from the selected one by more than 0,02 mA, the figure on the display will blink.

In case of the disconnection of the cable from the needle or the ECG electrode, (otherwise in case of the extraction of the needle from the patient skin,) the current will turn to zero, the display blinks and the stimulator generates a continuous sound.

If the needle has not been inserted in the patient yet, after four seconds the figure on the display will turn to zero (supplied current) and will start blinking but without the alarm sound.

On the upper left part of the display, the battery indicator will illuminate and blink if the battery charge starts to be insufficient to grant a proper functioning of the stimulator.

Since when the low battery indicator lights on the first time, the device can still operate for at least two hours of continuous utilization.

The display shows the selected parameters and which is the active selector for current / charge adjustment (knob or side control).

The "P" letter lights on when the stimulator automatically passes to "Precision" $\mu\text{s}/\text{MODE}$ (see par. 2.7.6).

2.6 Side Control (alternative to rotating knob for adjustments)

The side control buttons can be activated by pressing and holding on the arrow for at least 2 sec.

Its activation excludes the rotating knob from the possibility to adjust the current / charge: on the display this condition is highlighted by the disappearing of the symbol and the appearing

of the symbol.

The side control takes over the function of the rotating knob: each pressure on the arrow part of it corresponds to a clockwise step of the knob; each pressure on the lower part of it corresponds to an anticlockwise step of the knob. Keeping the lower or the upper button of the side control pressed, the current / charge will be adjusted continuously.

The utilization of the side control is particularly useful when the stimulator is used within a sterile field, with its protection bag (supplied separately), or for one hand manipulation.

Once the side control is activated, in order to deactivate it the stimulator must be switched off.

2.7 Safety features

2.7.1 Current self limitation at switching on

For safety reasons, the stimulator will not permit excessive currents before the needle has been inserted into the patient. Therefore, independently from the selected duration of the impulse, it will not be possible to initially set a current higher than 2,0 mA (or an equivalent charge).

Once the needle has been inserted into the patient, it's possible to select currents higher than 2,0 mA, as described in par. 2.4.

2.7.2 Sound tones

The stimulator generates a sound synchronized to the electrical impulse supplied to the patient (then at the selected frequency)

The tone of this sound remains unchanged for any current higher than 1,0 mA (or equivalent charge values)

When the current setting is reduced the sound increases by one tone each 0,1 mA of reduction up to 0,5 mA (highest tone).

For currents lower than 0,5 mA (or equivalent charge), the reached tone keeps unchanged.

From the menu of set up modification (point 2.3) it is possible to deactivate the sound during normal functioning. In that case the symbol will be crossed.

If the sound tone has been deactivated, check the display to make sure that the current is effectively supplied to the patient.

It is not possible to deactivate the alarm sounds.

2.7.3 Set / supplied current

The stimulator provides continuous and automatic check and adjustment of the supplied current to the selected one. The visualized current is then the selected one but also the supplied one, apart from the first four seconds from the variation of the selected current. (see par. 2.5).

In case of peaks of resistance preventing the stimulator from actually supplying the selected current, the visualized figure will blink* and will be the real one supplied.

If the needle has not been introduced into the patient yet, after four seconds the figure on the display will turn to zero (supplied current) and start blinking, but without any sound alarm.

2.7.4 Disconnection alarm

In case of disconnection **from** the needle or the ECG electrode, the displayed current turns to zero and blinks, and the stimulator generate a continuous tone.

2.7.5 “Safety” button

At any moment it's always possible to turn the current / charge immediately to zero by pressing the “Safety” button.

The pressure of the “Safety” button, further than turning to zero the supplied current, activates also a protection circuit turning the high tension generator circuit off and short circuits the output of the stimulator; to re-activate all circuits, the knob must be turned clockwise.

Pressing only the “Safety” button keeps all other set parameters unchanged.

2.7.6 “Precision” feature

The search for the nerve does not require an absolute precision. Therefore the stimulator qualifies higher variation steps (0,1 mA) for the adjustment of the current / charge (see par. 2.4) in order to avoid that too much turns of the knob are required for obtaining starting current values.

On the contrary, the fine adjustment of the needle tip position requires reduced values of the charge / current. Therefore the stimulator automatically qualifies a fine regulation (steps of 0,02 mA) if the current is set lower than 0,5 mA (or equivalent charge, as per table I par. 2.4).

When the “Precision” feature is activated, the letter P is shown on the display.

2.7.7 Exhausted battery

On the upper left part of the display, the battery indicator will illuminate and blink if the battery charge starts to be insufficient to grant the proper functioning of the stimulator.

However when the low battery indicator lights when the on button is activated, the device can still operate for at least two hours of continuous utilization.

At this time the battery should be changed as soon as possible, and before proceeding with further blocks.

2.7.8 Self switching off

If the set current is set at zero for 15 minutes without any further adjustment on any stimulator selector, the stimulator will automatically switch off.

3. INSTRUCTIONS FOR USE

3.1 Use of the device

The stimulator must be used by qualified personnel and only after completely reading the enclosed manual of instructions for use. These personnel must also have a sound knowledge of anatomy and neuro-physiology, and knowledge of the protocols in force for the execution of peripheral blocks.

The recommended approach is the following:

- start the search with impulses of relatively high intensity: for most cases an impulse of 300 μ s of duration with a charge of 300 ÷ 600 nC (1.0 ÷ 2.0 mA) is recommended;
- once the searched nerve has been identified, decrease progressively the intensity of the stimulus up to 90 - 150 nC (0.3 - 0.5 mA) in order to improve the position of the needle tip;
- if necessary, the precision of the positioning of the needle tip can be further improved by decreasing the current / charge (in this case the "Precision" feature will be automatically activated), or by selecting shorter impulse durations (100 or 50 μ S)
- once the required muscle contraction has been identified with the appropriated charge, following aspiration with a syringe, it's possible to inject a test dose of anaesthetic drug (around 1 ml).

The injection of anaesthetic drugs causes a variation of the nerve response to the stimulation; therefore after the injection it will not be possible to maintain the stimulation at the same levels of current.

Always test for aspiration of blood before any injection.

Using lower charge allows further precision but implies higher risks of direct contact with the nerve.

Once verified that the test dose has been successful, inject the full dose required.

Search with the same method all components of the nerve necessary for the complete execution of the block.

Since the device is not supposed to be sterile, do not bring it in a sterile field nor get it in touch with sterile personnel.

3.2 Check of the functioning

When the stimulator is switched on, all default parameters appear on the display:

300 μ s relevant to the selection of a duration of 300 μ s;

2 Hz relevant to the selection of a frequency of 2 Hz;

mA relevant to the selection of measure unit;

The figure shown on the LCD display will be 0000 if the selected unit is nC or 0.00 if the selected unit is mA.

If at the time of switching on the display is any different from what is described, the stimulator is not suitable for the use.

If the symbol of battery is blinking, this means that the battery is not charged enough and must be changed.

With the cable connected to the stimulator it is possible to verify that the device is working properly. To do so, adhere to the following procedure:

- to adjust the current or the charge; rotate the knob, verify that the value seen on the display is the desired value/setting. The screen will return to zero after about four seconds. By rotating the knob rapidly clockwise and anti-clockwise, you can linearly adjust the current setting; by rotating the knob slowly, a single click at a time, either clockwise or anti-clockwise, the corresponding value is displayed on the screen; verify that the sound tone which indicates current supply is clearly audible, and there is no flashing symbol of an exhausted battery;
- connect the output male plug of the cable (black wire) to the crocodile connection of the input cable (red wire); verify that, adjusting the current or charge by rotating the knob, the figure on the display is the selected one, after four seconds from the rotation of the knob, and stays still (not blinking).

If the result of the check has been positive, the stimulator can be used safely.

3.3 Conditions of use, cleaning and maintenance

THE PLEXYGON STIMULATOR DOESN'T NEED ANY TYPE OF PREVENTIVE MAINTENANCE.

In order to evaluate the proper functioning of the stimulator just follow the instructions in paragraph. 3.2. When the stimulator is not in use, it's recommended to switching it off, so as to preserve the battery charge. After use, the stimulator and its accessories have to be stored in the case

If the stimulator will not be used for a long period of time, it's recommended to remove the battery and store it also in the case. If the battery symbol is blinking, it means that the battery is exhausted and must be replaced as soon as possible; it's recommended to avoid the use of the stimulator in battery exhausted condition. The exhausted battery must be removed from the stimulator and discarded.

According to the electrical medical devices security standard (EN 606011/A2), we recommend to use the stimulator at a operating temperature between 10°C and 40°C and a relative humidity included between 30% and 75%.

The stimulator must be cleaned by a soft cloth dampen with natural water; it's possible to use also degreasing or sterilizing substances. The usage of ethylic alcohol is NOT allowed. Don't use sprays nor immerse the stimulator in any liquid.

A spilling or a stagnation of liquid or humidity on the device can cause a malfunction of the stimulator, therefore it is important to avoid any spills of liquid on the stimulator. In case of accidental spilling mop carefully any eventual residue.

Any technical modification of the device is not allowed, technical interventions on the device can be performed by Vygon Italy authorized personnel only to which the device must be sent in case of malfunctioning.

It's possible to sterilize the output cable by Ethylene Oxide sterilization at a temperature not higher than 50°C.

3.4 Accessories

The stimulator has been designed for performing peripheral nerve blocks. To this purpose, we have foreseen the utilization of the following accessories:

1. Electro-stimulation needle for the search of the nerve and the injection of the anaesthetic drugs. The stimulation needle must be connected to the stimulator's black cable. This needle must comply to the Essential Requirements foreseen by the 93/42/CEE Directive on the biocompatibility and is not supplied with the PLEXYGON stimulator. For safety and conformity reasons the approved needles are the "Locoplex" and "Echoplex" needles manufactured by Vygon S.A.- Ecoeu (France).
2. Continuous Peripheral Nerve Block: similarly to the previous point, similarly indwelling catheters must also comply with the same Essential Requirements. The approved sets are the "Techniplex", "Multiplex" and "Silverstim" manufactured by Vygon S.A. - Ecoeu (France).
3. Electrode ECG type, pre-gelled, to be connected to the red cable of the

stimulator. These electrodes are commonly available in the hospital, and must comply with the guide ISO 10993-1 on the compatibility, and are not supplied with the stimulator.


4. TECHNICAL DATA

4.1 General information


Power supply:	internal alkaline battery 9V, 500 mAh (6LR61);
Battery duration:	more than 25 hours of continuous functioning;
Impulse frequency:	1 Hz, 2Hz, 4 Hz, selectable by the “ μ S/MODE” button on the front panel.
Impulse duration:	50 μ s, 100 μ s, 300 μ s, selectable by the “MS/MODE” button on the front panel.
Current impulse:	adjustable by the knob on the front panel or the side control: from 0.0 mA to 6.0 mA for duration of 50 μ s; from 0.0 mA to 5.0 mA for duration of 100 μ s; from 0.0 mA to 4.0 mA for duration of 300 μ s.
Quantity of charge:	adjustable by the knob on the front panel or the side control: from 0 to 300 nC for duration of 50 μ s, from 0 to 500 nC for duration of 100 μ s, from 0 to 1200 nC for duration of 300 μ s.
Measure unit selection:	mA or nC, selectable by the menu.
Display:	LCD custom 4 digits and icons for the visualization of the current, or the quantity of charge, and for the set up of the stimulator.
Operating temperature:	Room temperature between 10°C and 40°C and relative humidity includes between 30% and 75%.

4.2 Output waveform

The output waveform has the typical course shown in the following figure, where two cases with resistance of 1000 Ω e 500 Ω are reported, with an impulse duration of 300 μ s and a current of 1,0 mA:

The maximum supplied current is the function of the resistance and of the duration selected, as shown in the following figure:

4.3 Precision

The total error on the current generation is of 0.2% + 5mV; it remains then always lower than 0.02 mA in the whole range of variation.

The error on the measurement of the supplied current is of around 0.01 mA in the whole range 0 - 6 mA. Therefore, when the user sets a current, the generated current has a precision better than 0.02 mA but, due to the following measurement error, the processor could even read a current with a difference of 0.02 mA as regards to the set current.

Then, the alarm of delivered current different from the set current, shown by the blinking value on the LCD display, starts for differences higher than 0.02 mA, in order to avoid fault readings.

4.4 Classification

Internally powered device whose applied part has been classified as Type BF. Classification according to IEC 60529 norm: IPX0, i.e. not protected against the ingress of liquids, degree of protection against the ingress of solid bodies not evaluated. Not suitable for the use in presence of anaesthetic mixtures flammable with air or oxygen or nitrogen protoxide, suitable for continuous functioning.

4.5 Liability

The manufacturer is responsible for electrical safety of the device in normal use, in usual conditions and in condition of first failure.

The liability on the safety, reliability and performances is only valid if:

- the stimulator has not been tampered in any way, and all adjustments, modifications or repairs have been carried out by authorized personnel only;
- the stimulator has been used in conformity of the present instructions for use.

4.6 Waste disposal and conformity to the WEEE directive

The following information is supplied pursuant to the LEGISLATIVE DECREE of 14 March 2014, no. 49, 'Implementation of the 2012/19/EU Directive on waste electrical and electronic equipment (WEEE)', concerning the reduction in the use of hazardous substances in electrical and electronic equipment as well as waste disposal.

The crossed-out refuse bin symbol, printed on the device or its package, indicates that at the end of its useful life the product must be collected separately from other waste.

The separate collection of this device at the end of its useful life is organised and managed by the manufacturer. The user who wants to dispose of this equipment will therefore need to contact the manufacturer and follow the system the latter has adopted to allow the separate collection of the device at the end of its useful life. An appropriate separate collection, leading to the decommissioned equipment's subsequent recycling, treatment and environmentally compatible disposal, will contribute to avoid possible negative effects on health and the environment, and promotes the re-use and/or recycling of the materials contained in the device.

Illegal disposal of the product by the user will result in the application of administrative sanctions pursuant to existing legislation.

Outside of the EU: follow local regulations for the disposal of medical and electronic equipment, including batteries.

4.7 Transport and storage

The device and its components do not need any particular temperature conditions. However, before using the instrument, check that the operating temperature corresponds to the point 3.3 of this manual.

4.8 Batch and date of manufacturing

The year of manufacturing is inserted in the batch number (see label on the battery room cover of the device).

The batch of manufacturing is composed by 6 (six) characters (two letters and 4 numbers), respectively and as follows subdivided:

- the first letter: the type of device (P stays for Plexygon)
- the second letter: the variant (A stays for: principal variant of Plexygon)
- the four numbers indicate:
 - the first two numbers: the year of the batch's manufacturing (ex. 10 stays for year 2010)
 - the third and fourth number: week of manufacturing (ex. 06 stays for first week of February).

4.9 Electromagnetic compatibility

Table I - Manufacturer's declaration - Electromagnetic emissions

The PLEXYGON is intended for use in the electromagnetic environment specified below. The customer and/or the user of the PLEXYGON should assure that it is used in an electromagnetic environment as described below.

Emissions test	Compliance	Electromagnetic environment guidelines
RF emission CISPR 11	Group I	The PLEXYGON uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The PLEXYGON is suitable for use in all establishments, including domestic establishments and those directly connected to the public power supply network that supplies buildings used for domestic purposes.

Table 2 - Manufacturer's declaration - Electromagnetic immunity - Enclosure port

Phenomenon	Basic EMC standard or test method	Immunity test levels required for Professional healthcare facility environment	Compliance test levels
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Radiated RF electromagnetic fields	IEC 61000-4-3	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See Table 4	See Table 4
Rated power frequency magnetic fields	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz

Table 3 - Manufacturer's declaration - Electromagnetic immunity - Patient coupling port

Phenomenon	Basic EMC standard or test method	Immunity test levels required for Professional healthcare facility environment	Compliance test levels
Electrical fast transients / bursts	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V 0,15 MHz – 80 MHz 6V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3V 0,15 MHz – 80 MHz 6V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz

Table 4 - Test specifications for enclosure port immunity to RF wireless communications equipment according to IEC 61000-4-3

Test Frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity Test level (V/m)
385	380 - 390	TETRA 400	Pulse modulation 18 Hz	1,8	0,3	27
450	430 - 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0,3	28
710	704 - 787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0,3	9
745						
780						
810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0,3	28
870						
930						
1720	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0,3	28
1845						
1970						
2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0,3	28
5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0,3	9
5500						
5785						

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the PLEXYGON, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

5. REGULATORY ISSUES

5.1 Applied standards

The following are the International regulations, European Union Directives and International Technical Standards that were followed to design, build and test Plexygon.

- Directive 93/42/EEC and following amendments (2007/47/EEC included)
- EN ISO 13485 Medical equipment - Quality system management
- EN ISO 15223-1 Symbols to be used with medical device labels, labelling and information to be supplied
- EN ISO 14971 - Medical devices - Application of risk management
- MED.DEV 2.12/1 Guidelines on a medical device vigilance system
- MED.DEV 2.7/1 Guidelines on Clinical Evaluation
- IEC 60601-1:2005 + Ec1:2006; IEC 60601:2006 + EC:2010 Medical electrical equipment – General safety standards
- IEC 60601-1-2 and EN 60601-1-2 Medical electrical equipment part 1: general requirements for safety. Collateral standard: electromagnetic compatibility - requirements and tests
- EN 60601-1-4 Medical electrical equipment part 1-4: general requirements for collateral standard - Programmable electrical medical systems.
- EN 60601-2-10/A1 Medical electrical equipment part 2: particular requirements for the safety of nerve

5.2 Declaration of conformity

VYAGON ITALY S.R.L. declares that the device PLEXYGON complies with 93/42/EEC Directive and further amendments included in 2007/47/EEC concerning medical devices, according with the Essential Requirements of the Annex I. Therefore the medical device PLEXYGON has the CE marking CE 0123 (see below) according with Annex XII of the 93/42/EEC Directive.

6. WARRANTY

Vygon Italia s.r.l. guarantees that, for all failures deriving from manufacturing defect during the period of the warranty, the stimulator will be repaired (workmanship and spare parts) free of charge for the customer.

For the period of the warranty please refer to the local Vygon subsidiary

Warranty does not apply in case of failures due to falls, crashes, negligence, unsuitable or unauthorized maintenance, or anyhow attributable to a misuse of the device.

Warranty is excluded also for wear and tear parts as batteries, cables, etc. Maintenance, under warranty or not, should be carried out exclusively by Vygon Italia s.r.l. or authorized centres.

Warranty is excluded in case of interventions carried out by unauthorized thirds.

Vygon Italia or the authorized center reserve the right to decide if to repair or replace the defective parts or stimulators. In case of replacement, the replaced parts will remain of property of Vygon Italia.

In order to require warranty repairs, a sale invoice or other relevant document and the valid serial number of each device must be presented.

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- Directive 93/42/EEC and following amendments (2007/47/EEC included)
- EN ISO 13485 Medical equipment - Quality system management
- EN ISO 15223-1 Symbols to be used with medical device labels, labelling and information to be supplied
- EN ISO 14971 - Medical devices - Application of risk management
- MED.DEV 2.12/1 Guidelines on a medical device vigilance system
- MED.DEV 2.7/1 Guidelines on Clinical Evaluation
- IEC 60601-1:2005 + Ed1:2006; IEC 60601:2006 + EC:2010 Medical electrical equipment – General safety standards
- IEC 60601-1-2 and EN 60601-1-2 Medical electrical equipment part 1: general requirements for safety. Collateral standard: electromagnetic compatibility - requirements and tests
- EN 60601-1-4 Medical electrical equipment part 1-4: general requirements for collateral standard - Programmable electrical medical systems.
- EN 60601-2-10/A1 Medical electrical equipment part 2: particular requirements for the safety of nerve

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- Directiva 93/42/EEC y posterior modificación (2007/47/CE incluida)
- EN ISO 13485 Medical equipment - Quality system management
- EN ISO 15223 Graphical symbols for use in the labelling of medical devices
- UNI EN ISO 14971 - Medical devices - Application of risk management
- MED.DEV 2.12/1 Guidelines on a medical device vigilance system
- MED.DEV 2.7/1 "Guidelines on Clinical Evaluation"
- EC 60601-1:2005 + Ec1:2006; IEC 60601:2006 + EC:2010 Medical electrical equipment – General safety standards
- IEC 60601-1-2 ed EN 60601-1-2 Medical electrical equipment part 1: general requirements for safety. Collateral standard: electromagnetic compatibility - requirements and tests
- EN 60601-1-4 Medical electrical equipment part 1-4: general requirements for collateral standard - Programmable electrical medical systems.
- EN 60601-2-10/A1 Medical electrical equipment part 2: particular requirements for the safety of nerve

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VYGON ITALIA Srl ▼ GRUPPO VYGON