REGIONAL ANESTHESIA
Peripheral blocks

User manual / Mode d'emploi / Manuale d'uso
Instrucciones de uso / Gebrauchsanweisung / Manual do usuário
Brugervejledning / Bruksanvisning / Upute za upotrebu
Instrucții de utilizare / Manual de instruções

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Contents

1 GENERAL	2
1.1 Intended use	2
1.2 Indications	2
1.3 Warnings and precautions	3
1.4 Contraindications	4
1.5 Undesirable effects	4
2 OPERATION OF THE STIMULATOR	4
2.1 Preparing the device for use	4
2.1.1 Checking the parts list	4
2.1.2 Installation/replacement of battery	4
2.1.3 Cable connection	4
2.2 Switching the device on/off	5
2.3 Changing operating parameters	5
2.4 Selecting charge or current	5
2.5 Display	5
2.6 Side button	6
2.7 Safety functions	6
2.7.1 Self-limiting current at power-on	6
2.7.2 Sound tone	6
2.7.3 Current set/ delivered	6
2.7.4 Disconnection warning sound	6
2.7.5 "Safety" button	6
2.7.6 "Precision" Function	6
2.7.7 Low battery	7
2.7.8 Auto shutdown	7
3 METHOD OF USE	7
3.1 Using the device	7
3.2 Verification of operation and safety before each use	7
3.3 Conditions of use, maintenance and cleaning	7
3.4 Medical devices for use in combination with PLEXYGON	8
4 TECHNICAL DESCRIPTION	9
4.1 General Information	9
4.2 Output waveform	9
4.3 Accuracy	9
4.4 Classification	10
4.5 Liability	10
4.6 Disposal and compliance with the WEEE Directive	10
4.7 Environmental conditions for transportation and storage of the device	10
4.8 Electromagnetic compatibility	11
4.9 Essential performance	12
5 SERIOUS INCIDENT REPORTING	12
6 WARRANTY	13
7 DESCRIPTION OF SYMBOLS ON THE LABEL	13

1 GENERAL

The PLEXYGON medical device has been manufactured according to current international safety standards. Please read this manual and all warnings carefully before use, as recommended by the symbol on the device label.

PLEXYGON's IFUs are also available in electronic format at www.vygon.it.

NOTE: The use of the stimulation needles/catheters required for stimulation is described in the accompanying instructions for use, which are not part of this user manual.

1.1 Intended use

PLEXYGON is an electrostimulator for localisation of peripheral plexuses and nerves in regional anaesthesia procedures or pain therapy.

Nerve localisation is achieved by generating appropriate electrical stimulation, through a dedicated needle electrode, then inducing a muscle reaction on the afferents of the nerve reached by the stimulation. Identifying a specific nerve allows the anaesthetic to be injected in close proximity to it, enabling local anaesthesia before surgery or for pain therapy. The device is intended for use in hospitals, laboratories, clinics and medical practices.

INTENDED PATIENTS

Peripheral nerve block with electrical stimulation is typically indicated for patients requiring limb surgery: peripheral nerve block guided by electrical stimulation can be used to temporarily block the sensation of pain in a specific limb, or part of a limb, during surgery; it is therefore often used in orthopaedic operations, such as hand and lower limb procedures.

Furthermore, electrical stimulation-guided anaesthesia is a technique mainly used for patients who cannot or do not want to undergo general anaesthesia, or for whom traditional anaesthetic drugs are contraindicated. The use of electrical stimulation for peripheral regional anaesthesia is indicated for all patients for whom this type of anaesthesia is considered clinically preferable to general anaesthesia, except for breastfeeding women and patients with pacemakers or other implanted or connected stimulators, in whom specialist advice is required. For paediatric patients, ultrasound guidance or a combined technique (ultrasound + electrical stimulation) is generally preferred.

Not recommended for use in pregnant women and newborns.

INTENDED USERS

The intended user is an anaesthesiologist experienced in peripheral plexus and nerve blocks, with the necessary knowledge of anatomy and neurophysiology, as well as the relevant protocols for performing such blocks.

EXPECTED INDIRECT CLINICAL BENEFITS

- **Precise guidance for needle insertion:** the electrostimulator helps precisely locate the target nerve. Electrical stimulation produces muscle contractions in the target area for regional anaesthesia, confirming the proximity of the needle to the relevant nerve and thereby improving the effectiveness of the block.
- **Effective pain management:** electrical stimulation-guided peripheral nerve blocks provide better nerve identification allowing more targeted analgesia, reducing intraoperative and postoperative pain, and improving patient comfort
- **Reduction in anaesthetic use:** by allowing precise identification of the nerve location and enabling injection in close proximity to it, electrical stimulation makes it possible to limit the amount of anaesthetic required, thereby reducing the risk of overdose.
- **Reduced complications:** the use of the electrostimulator decreases the risk of incorrect needle insertions, minimising the risk of nerve or vascular injury and other complications related to the blocking technique.
- **Faster postoperative recovery:** peripheral nerve blocks (PNBs) have been associated with fewer systemic side effects and improved pain control, contributing to a faster return to normal activities.

1.2 Indications

The device is indicated for identifying a specific nerve or nerve bundle and for enabling injection of anaesthetic or analgesic agents in close proximity to the target area, thereby facilitating regional anaesthesia before surgery or for pain management.

For greater safety and efficacy in locating the nerve, it is preferable to use nerve stimulation in combination with ultrasound guidance. The combined use of the two techniques allows direct visualisation of the anatomical structures and functional confirmation of the needle position, reducing the risk of neurological complications and improving the accuracy of the block. Combined use is particularly recommended in the paediatric population.

1.3 Warnings and precautions

- If PLEXYGON is not working as described in this manual, do not use the device and seek technical assistance.
- The combined use (dual guidance) of ultrasound and electrical stimulation in peripheral nerve blocks with a motor component in adult patients increases the efficacy and safety of the surgical block, reducing the risk of serious neurological complications.

- The user must check the safety, status of the device and check for damage or malfunction before use. Carry out verifications as specified in Section 3.2.
- If, when switched on, the display differs from what is described in Section 3.2, the stimulator must not be used.
- If the battery symbol flashes, the battery is low and must be replaced before starting the nerve search procedure.
- The injection of anaesthetic or analgesic agents alters the nerve response to stimulation; therefore, after injection it will not be possible to maintain stimulation at the same current levels.
- Using a lower charge allows further precision but implies higher risks of direct contact with the nerve.
- If the acoustic signal is disabled, the display should be checked to verify that current is actually being delivered to the patient (section 2.7.2).
- PLEXYGON is not sterile; do not place it in a sterile field or allow it to come into contact with sterile personnel without appropriate precautions.
- The use of the stimulator can interfere with the functioning of pacemakers or other implanted or connected stimulators. Therefore do not use the stimulator without consulting a specialist in these devices.
- Do not use the stimulator in the presence of flammable gases.
- Do not use the stimulator in case of malfunction, leakage of battery fluids or accidental spills onto the stimulator.
- Do not use non-original medical devices intended for use in combination with PLEXYGON. The use of non-original medical devices intended for use in combination with PLEXYGON may alter the characteristics of the electrical signal, increase electromagnetic emissions and/or reduce the electromagnetic immunity of the device, thereby compromising its safety features. Check the compatibility between PLEXYGON and other medical devices intended for use in combination with it. Handle with particular care (Section 3.4).
- Do not use near microwave and/or short wave therapy devices.
- High risk of cardiac fibrillation due to application of electrodes on the chest.
- The stimulation should not be applied on or through the head, directly on the eyes, covering the mouth, on the front of the neck (especially the carotid sinus), or from electrodes placed on the chest, the upper back or over the heart.
- Do not touch the battery contacts and the patient at the same time.
- To reduce the risk of battery leaks, remove the battery from the device when not in use for a long time.
- Spillage or accumulation of liquid or moisture on the device may cause the stimulator to malfunction; therefore, it is important to avoid any liquid spills on the stimulator. In case of accidental spillage, dry the area immediately.
- To achieve optimal nerve stimulation, always ensure the electrode used is intact and not dry.
- To avoid damaging the connecting cables and the device, do not hold or carry the device by its connecting cables and/or other medical devices intended for use in combination with PLEXYGON. Do not wrap the cable around the device or other equipment.
- For safety reasons, never operate PLEXYGON if the battery leaks. Return the device to Vygon Italia authorised personnel for proper cleaning. The same applies in the event of liquid leakage.
- Do not replace the battery if the patient is connected to the device or the device is turned on.
- If the device is not in use, switch it off and remove the battery.
- For device disposal please refer to the instructions specified in Section 4.6 of the user manual.
- No modifications to the PLEXYGON device are permitted; technical interventions may only be carried out by Vygon Italia personnel.
- Electrodes exposed to current densities higher than $2\text{mA}/\text{cm}^2$ may require special attention from the operator. Refer to the IFU supplied with the electrodes.

WARNING: The simultaneous connection of the patient to a surgical high-frequency device may cause burns at the stimulator electrodes or damage to the stimulator itself.

WARNING: Mobile, radio frequency communication devices (including cables for antennas and external antennas) must not be used closer than 30 cm (12 inches) to any part of the device, including the cables specified by the manufacturer. Otherwise, the performance of the device may be affected.

WARNING: Verify that the cables do not come into contact with other cables or directly with the patient.

WARNING: Use of this device near or over (above or below) other devices is prohibited as it may cause improper functioning. If this type of use is necessary, the device in question and the others involved must be closely monitored to verify their correct functioning.

WARNING: Pay attention to the position of metallic implants in the tissue (e.g. plates or electrode leads), which may channel stimulation signals to other sites and cause harmful effects. Implanted electronic equipment can be damaged by the stimulation current, which in turn may lead to malfunction of the implants or even their destruction.

WARNING: The cable supplied may only be connected to PLEXYGON and to other medical devices intended for use with it (stimulation needle, continuous peripheral block catheter, and ECG-type electrode applied to the patient).

1.4 Contraindications

The use of PLEXYGON is contraindicated in the following cases:

- Active implantable device: use of the stimulator may interfere with the functioning of pacemakers or other implanted or connected stimulators. Therefore, do not use the stimulator without consulting a specialist in these devices.
- Non-cooperative patients: electrical stimulation requires the patient's active cooperation, as the anaesthetist may request feedback on sensation or on the extent of the anaesthetic block. Patients who are unable to communicate or cooperate adequately may not be suitable candidates.
- Pregnant women and newborns
- Pre-existing neurological disorders: in some pre-existing neurological conditions, such as polyneuropathy or multiple sclerosis, electrical stimulation may be inadvisable or require further evaluation. These conditions may influence sensitivity to stimulation and the results of the procedure.

It is essential that the anaesthetist carefully assesses the patient's clinical situation and medical history before administering peripheral regional anaesthesia, irrespective of whether electrical stimulation is used or not (e.g. local infection, allergy to anaesthetics, bleeding disorders, skin lesions or defects, etc.). Only a qualified physician can provide an adequate assessment of individual risks and benefits.

1.5 Undesirable effects

Electrostimulators are used in peripheral regional anaesthesia to identify and block specific nerves, but like all medical procedures, they can have undesirable effects, such as:

- Pain, discomfort, or paraesthesia: During needle insertion and electrical stimulation, the patient may experience pain, discomfort, or abnormal sensations, such as tingling or burning.
- Nerve injury: needle insertion and inappropriate use of electrical stimulation can cause nerve injury, which may lead to temporary or, in rare cases, permanent neuropathy.
- Nausea and/or vomiting
- Vascular perforation: possible if the needle is inserted close to vascular structures.
- Pneumothorax and Horner syndrome: may occur following blocks performed near the thoracic cavity or cervical plexus.
- Local anaesthetic systemic toxicity (LAST)

These side effects vary depending on the technique used, the physician's experience, and the patient's individual response. It is important that the procedure be performed by an experienced anaesthesiologist and in a controlled environment to minimise risks.

2 OPERATION OF THE STIMULATOR

2.1 Preparing the device for use

2.1.1 Checking the parts list

Verify that the case contains the following components in addition to the device:

- 1 cable for connecting the stimulator to the needle and electrode on the patient
- 1 9V alkaline battery - 500 mAh (type 6LR61)
- IFU of the device

2.1.2 Installation/replacement of battery

Open the cover of the battery compartment, located on the lower rear of the device, by gently pressing the latch and lifting it.

At this point, connect the 9V-500mAh alkaline battery (type 6LR61) to the appropriate clip inside the compartment, ensuring correct alignment of the positive and negative terminals, ensuring correct alignment of the positive and negative terminals and checking that the clip is securely in place. Then insert the battery into the compartment and close the cover by following the opening procedure in reverse order.

2.1.3 Cable connection

Once the battery compartment is closed, connect the stimulator to the needle and to the ECG-type electrode placed on the patient. Insert the three-pin plug of the cable into the three-pin socket on the upper side of the device, then connect the positive terminal (red alligator clip) to the ECG-type electrode and the negative terminal (black socket) to the needle.

2.2 Switching the device on/off

To switch on the device, press the on/off button located on the front panel. When the device is switched on, the display shows the default settings: the current (0.00 mA at start-up) associated with each generated pulse, the unit of measurement (mA), the frequency (2 Hz) and the stimulus duration (300 μ S).

Before proceeding further, it is recommended check that the device is operating properly, as described in Section 3.2 below.

To switch off the device, press and hold the button for 2 seconds.

2.3 Changing operating parameters

To change the operating parameters, press and hold the μ S/MODE button at least 2 seconds. Once the menu is open, the editable parameter is highlighted; turn the knob to change its value. Pressing the μ S/MODE button again selects the next parameter. The selection sequence is as follows:

Frequency:	Selectable options:	2 - 1 - 4 Hz (default 2 Hz)
Duration:	Selectable options:	300 - 100 - 50 μ S (default 300 μ S)
Units of measurement:	Selectable options:	mA – nC (default mA)
Sound:	Selectable options:	on – off (default: on)

A further press of the μ S/MODE button exits the parameter selection menu and returns the device to normal operation.

During operation, a short press and immediate release of the μ S/MODE button allows the pulse duration to be changed while keeping the delivered charge (expressed in nC) constant. When displayed in mA, the current will vary to maintain a constant delivered charge (e.g. duration 300 μ S - current 0.5 mA; duration changed to 100 μ S - current 1.5 mA; duration 50 μ S - current 3.0 mA).

2.4 Selecting charge or current

The knob on the upper right of the front panel allows adjustment of the current or charge delivered by the device with each pulse.

With each knob rotation, the set current changes by 0.1 mA or 0.02 mA, depending on whether the set current is greater or less than 0.5 mA, respectively. Turning the knob clockwise increases the set current, while turning it counterclockwise decreases it. During the fine adjustment phase, the symbol "P" appears on the display.

NOTE: When the display unit is "nC," the increment in charge value corresponds to the product of the current and the pulse duration, according to the values shown in the following table:

Table 1: Charge variation step Q (nC) for different pulse durations


	300 μ S	100 μ S	50 μ S
Q $\geq$ 150 nC	30 nC	-	-
Q < 150 nC	6 nC	-	-
Q $\geq$ 75 nC	-	10 nC	-
Q < 75 nC	-	2 nC	-
Q $\geq$ 25 nC	-	-	5 nC
Q < 25 nC	-	-	1 nC

The adjustment range of the set current by turning the knob is as follows: 0 to 4 mA for pulse durations of 300 μ s (i.e. 0 to 1200 nC), 0 to 5 mA for pulse durations of 100 μ s (i.e. 0 to 500 nC) and 0 to 6 mA for pulse durations of 50 μ s (i.e. 0 to 300 nC).

During operation, each pulse is accompanied by a sound. This sound increases in pitch every 0.1 mA between 0.5 and 1.0 mA.

In any case, under no circumstances and for no pulse duration may a current value greater than 2.0 mA (or the equivalent charge) be set before the needle is inserted into the patient (see Section 2.7.1).

2.5 Display

The display is used both to show the current or charge associated with each pulse and to show selected parameters (pulse duration, frequency, display unit, selection button, battery status). The displayed parameter corresponds to the value actually delivered by the device, except when the adjustment knob is operated; in that case, the displayed parameter corresponds to the set value, both during adjustment and for the next four seconds.

If the delivered current differs from the set current by more than 0.02 mA, the value shown on the LCD display will flash.

If the needle or ECG electrode is disconnected, or if the needle is removed from the patient's skin, the current will drop to zero and the stimulator will emit a continuous sound. If the needle has not yet been inserted into the patient, after 4 seconds the

display will return to zero (delivered current), but without a warning sound.

A flashing battery symbol appears at the upper left of the display if the battery charge is insufficient for proper device operation.

From the time the low battery symbol first appears, the device remains capable of functioning for at least two hours of continuous use.

The display highlights the set parameters and indicates which button is active for current/charge adjustment.

The letter "P" lights up when the stimulator automatically switches to the "Precision" function (see Section 2.7.6).

2.6 Side button

The side button can be activated by holding down the top of the button for at least 2 seconds.

Activation of the side button deactivates the upper knob: this is indicated on the display by the disappearance of the symbol ☉ and the appearance of the symbol ◀▶.

The side button replaces the upper knob: pressing the upper button corresponds to a clockwise rotation of the knob; pressing the lower button corresponds to a counterclockwise rotation. Holding down the upper or lower buttons results in a continuous increase or decrease in current, respectively.

Use of the side button is particularly useful when operating the stimulator in a sterile field, with the appropriate protection sleeve (supplied separately), or for one-handed operation.

Once the side button is activated, the stimulator must be switched off to deactivate it.


2.7 Safety functions

2.7.1 Self-limiting current at power-on

For safety reasons, the device does not allow excessive currents to be set before the needle is inserted into the patient. Therefore, regardless of the pulse duration selected, a current greater than 2.0 mA (or equivalent charge) cannot be set at start-up.

Once the needle is inserted into the patient for the first time, higher currents may be selected, as outlined in Section 2.4.

2.7.2 Sound tone

The stimulator emits a sound pulse synchronous with the electrical pulse delivered to the patient (i.e. at the same frequency as the one set). The pitch of this sound remains constant for set currents above 1.0 mA (or equivalent values in terms of charge).

Below this level, the pitch increases by one tone for every 0.1 mA reduction in current, up to of 0.5 mA.

For currents below this value, the sound remains constant.

From the parameter adjustment menu (section 2.3) it is possible to disable the sound in normal operation; in this case the symbol will be crossed out.

However, it is not possible to switch off the warning sounds.

2.7.3 Current set/ delivered

The stimulator automatically checks and adjusts the delivered current with respect to the set current. Therefore, the displayed current is the set current and also the delivered current, except for the first 4 seconds after the change in the set current (see Section 2.5).

In case of resistance spikes that prevent the stimulator from delivering the set current, the value shown on the display will flash, corresponding to the actual delivered current.

If the needle has not yet been inserted into the patient, after 4 seconds the display will return to zero (delivered current), but without a warning sound.

2.7.4 Disconnection warning sound

If the needle or ECG electrode is disconnected, the current drops to zero, and the stimulator emits a continuous sound.

2.7.5 "Safety" button

At any time, it is possible to reset the current output or charge by pressing the "Safety" button. Pressing the "Safety" button, in addition to resetting the set current to zero, enables a protection circuit that turns off the high-voltage generating circuit and short-circuits the device output; restoration of the high-voltage generating circuit and of the output are achieved by turning the current control knob clockwise.

Pressing the "Safety" button does not change other parameter settings.

2.7.6 "Precision" Function

The nerve search phase does not require absolute precision, so the stimulator enables larger increments (0.1 mA) for current/charge adjustment (see Section 2.4) to avoid excessive knob rotations when reaching the required starting values.

Conversely, fine adjustment of the needle tip position requires small current/charge values; therefore, the stimulator automatically enables fine adjustment of the current (0.02 mA increments) for values below 0.5 mA (or equivalent charge values, as per Table 1 of Section 2.4).

Activation of the "Accuracy" function is indicated by the appearance of the letter "P" on the display.

2.7.7 Low battery

A flashing battery symbol appears at the upper left of the display if the battery charge is insufficient for proper device operation.

From the time the low battery symbol first appears, the device remains capable of functioning for at least two hours of continuous use.

However, replace the battery as soon as possible, and in any case before performing blocks on other patients.

2.7.8 Auto shutdown

If the set current remains at zero without further button presses for 15 minutes, the stimulator will automatically switch off.

3 METHOD OF USE

3.1 Using the device

Before each use, check the functionality and safety of the device, as described in Section 3.2.

The display must be constantly monitored during operation to confirm the actual current delivery.

The recommended approach is as follows:

- Begin the search with relatively high-intensity pulses: in most cases, a stimulus of 300 μ s with a charge of 300 + 600 nC (1.0 $\pm$ 2.0 mA) is appropriate;
- Once the target nerve is identified, gradually reduce the stimulus intensity to about 90-150 nC (0.3-0.5 mA) to improve needle placement;
- If necessary, the accuracy of needle placement can be improved by further reducing the delivered current/charge values (in this case, the "Accuracy" function is automatically activated), or by switching to shorter pulse durations (100 or 50 μ s);
- Once the desired contraction has been identified with the appropriate amount of charge, a test dose of anaesthetic/analgesic (about 1 ml) can be injected.

Once the effectiveness of the test dose has been confirmed, inject the desired full dose.

Use this method to identify all nerve components required for complete block execution.

3.2 Verification of operation and safety before each use

Before each use, check the functionality of the display, beeper, knob, cable, and output pulse generation. For this test, proceed as follows:

- Connect PLEXYGON's cable to its front three-pin socket;
- Turn on the device and verify that for 1.5 seconds all segments are illuminated and the warning sound is audible;
- Adjust the current or charge by turning the knob, verifying that the value shown on the display is that set during the adjustment and remains for the next four seconds, then returns to 0 after four seconds. Verify that by turning the knob quickly, clockwise or counterclockwise, the desired current setting can be adjusted linearly; verify that by turning the knob slowly, one click at a time, clockwise or counterclockwise, the current values shown on the display correspond to the expected values. Verify that the sound tone signalling current delivery is clearly audible. Verify that the flashing low battery symbol is not present.
- Connect the male connector of the connection cable to the needle (black) to the alligator clip of the return cable (red); verify that by adjusting the current or charge with the knob, the value shown on the display is as set, even after the four seconds, i.e. the current generated by the device is the same as the set current (non-flashing display). Then disconnect the connection cable from the alligator clip of the return cable and verify that the current value returns to 0 and the stimulator emits a continuous sound.

If the test is successful, the stimulator can be used safely.

If the test is unsuccessful, the stimulator must not be used. Contact Vygon technical support to assess whether the cable needs to be replaced or the device returned for service.

3.3 Conditions of use, maintenance and cleaning

PLEXYGON does not require any kind of preventive maintenance; correct operation is verified by the user before each procedure using the test described in Section 3.2.

When the device is not in use, switch it off to preserve battery life.

After use, the device and its components should be stored in the case.

For long periods of non-use, remove the battery and store it in the compartment inside the case.

If, with the device on, the battery symbol at the upper left side of the display flashes, the battery is low and should be replaced as soon as possible. The device should not be used with a low battery. The discharged battery should in any case be removed from the device.

Referring to the safety standards for electromedical equipment (EN 60601-1/A2), it is recommended to use the stimulator at room temperature between 10°C and 40°C and at a relative humidity between 30% and 75%.

PLEXYGON and its connecting cable must be cleaned before and after each use. Cleaning should be carried out with a soft cloth moistened with water; degreasing or disinfecting substances such as glutaraldehyde may be used, but ethyl alcohol is not recommended. Do not use sprays or soak the stimulator.

3.4 Medical devices for use in combination with PLEXYGON

The stimulator is designed for peripheral nerve blocks. For this purpose, the following medical devices, not included in the case, are intended for use with PLEXYGON:

1. Electrical stimulation needle for nerve localisation and anaesthetic injection. The stimulation needle must be connected to the stimulator's black cable. This needle must comply with the Essential Requirements of Directive 93/42/EEC or Regulation (EU) 2017/745 on biocompatibility and is not supplied with the PLEXYGON stimulator. The needles validated by the manufacturer of PLEXYGON are Locoplex and Echoplex + produced by Vygon S.A. (France) and Vygoplex produced by MDL S.r.l. (Italy).
2. Stimulable continuous peripheral block catheter: As with stimulation needles, indwelling catheters must also comply with the same Essential Requirements. The approved sets are the Techniplex, Multiplex and Silverstim manufactured by Vygon S.A. (France).
3. ECG-type electrode, to be connected to the red cable of the stimulator. These ECG electrodes are commonly available on the EU market as medical devices and are not supplied with the stimulator.

4 TECHNICAL DESCRIPTION



4.1 General Information

DIMENSIONS	• max. width: 94 mm - Min. width: 58 mm • max. height: 39.5 mm - Min. height: 23 mm
WEIGHT	• Weight: 186 g • Weight with the battery: 232 g
POWER SUPPLY	Internal alkaline battery 9V, 500 mAh (6LR61)
NUMBER OF USES	Unlimited
STERILITY	Non sterile
SOFTWARE VERSION	4.00
BATTERY DURATION	25 hours of continuous functioning
CURRENT ACCURACY	$\pm 5\%$
IMPULSE FREQUENCY	1 Hz, 2 Hz, 4 Hz selectable with "µS/MODE" button on front panel
IMPULSE DURATION	50 µs, 100 µs, 300 µs selectable with "µS/MODE" button on front panel
CURRENT IMPULSE	• from 0.0 mA to 6.0 mA for duration of 50 µs • from 0.0 mA to 5.0 mA for duration of 100 µs • from 0.0 mA to 4.0 mA for duration of 300 µs
CHARGE	• from 0 to 300 nC for duration of 50 µs • from 0 to 500 nC for duration of 100 µs • from 0 to 1200 nC for duration of 300 µs
UNIT OF MEASUREMENT	mA or nC
DISPLAY	custom
OPERATING TEMPERATURE	Room temperature between 10°C and 40°C Relative humidity between 30% and 75%. Atmospheric pressure range from 70.0 to 106.0 kPa
MATERIALS IN CONTACT WITH THE USER	ABS
PROTECTION CLASS	internally powered
DEVICE TYPE	BF
APPLIED PARTS	positive terminal (red alligator clip) and negative terminal (black socket) The maximum surface temperature of the applied parts is 42°C.
IP RATING	IPX0
ABSORPTION	Less than 20 mA

4.2 Output waveform

The output waveform has the typical pattern shown in the following figure, where two cases of 1000 and 500 $\square$ load resistance are shown, with a pulse duration of 300 µs and a current of 1.0 mA.

The maximum current amplitude delivered is a function of the load resistance and pulse duration selected, as shown in the following figures:

4.3 Accuracy

The overall error on current generation is $\pm 5\%$, in the 0 - 6 mA range it averages 0.01 mA.

If the set current differs from the measured current, this is indicated by the value flashing on the display, when the deviation is at least 0.02 mA.

4.4 Classification

Internally powered device with applied parts classified as Type BF. Classification according to IEC 60529: IPX0, i.e. not protected against liquid ingress, degree of protection against solid ingress not evaluated. Not suitable for use in the presence of flammable anaesthetic mixture with air or oxygen or nitrous oxide, suitable for continuous operation.

4.5 Liability

The manufacturer is responsible for the electrical safety of the device in normal use, under usual conditions, and under first-fault conditions.

Liability for the purposes of device safety, reliability and performance applies only if:

- The device has not been tampered with in any way, and all adjustments, modifications, or repairs have been made by authorised personnel;
- The device is used in accordance with the operating instructions.

4.6 Disposal and compliance with the WEEE Directive

The following information is provided in accordance with LEGISLATIVE DECREE No. 49 of 14 March 2014 "Implementation of Directive 2012/19/EU on waste electrical and electronic equipment (WEEE)" regarding the reduction of hazardous substance use in electrical and electronic equipment and waste disposal.

The crossed-out wheeled bin symbol on the device and its packaging indicates that, at the end of its service life, the product must be collected separately from other waste.

Separate collection of this end-of-life equipment is organised and managed by the manufacturer.

Users who wish to dispose of the device must clean it in accordance with Section 3.3 and contact the manufacturer to follow the system in place for the separate collection of end-of-life devices.

Adequate separate collection of the decommissioned device for subsequent environmentally sound recycling, treatment and disposal helps to avoid possible adverse effects on the environment and health, and promotes the reuse and/or recycling of the materials from which the device is made.

Unauthorised disposal of the product by the user will result in administrative sanctions as provided for in current regulations.

Outside the European Community: follow local regulations for disposal of medical and electronic equipment, including batteries.

4.7 Environmental conditions for transportation and storage of the device

The device and its critical components do not require special environmental conditions beyond those of normal living environments.

However, before using the device, check that it meets the environmental conditions for use indicated in Section 3.3.

4.8 Electromagnetic compatibility

Table 1 - Manufacturer's Declaration - Electromagnetic Emissions

The PLEXYGON electrostimulator is intended for use in the electromagnetic environment specified below. The customer or the user of the PLEXYGON electrostimulator must ensure that it is used in such an environment.

Emission testing	Compliance level	Electromagnetic environment - Guidance
RF emissions CISPR 11	Group 1	The PLEXYGON electrostimulator uses RF energy only for its internal operation. As a result, its RF emissions are very low and not likely to cause any interference with electronic devices placed nearby.
RF emissions CISPR 11	Class B	The PLEXYGON electrostimulator is suitable for use in all environments, including domestic settings and those directly connected to a low-voltage public power supply intended for household use.

Table 2 - Manufacturer's declaration - Electromagnetic immunity - Enclosure port

Phenomenon/ Immunity test	Basic EMC Standards or test method	Level of immunity test required for professional environment (healthcare facilities)	Compliance level
Electrostatic discharges	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV in air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV in air
Radiated RF electromagnetic fields	IEC 61000-4-3	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz
Proximity fields from RF wireless communication devices	IEC 61000-4-3	See Table 4	See Table 4
Magnetic fields at grid frequency	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz

Table 3 - Manufacturer's declaration - Electromagnetic immunity - Patient cable port/Patient connection

Phenomenon/ Immunity test	Basic EMC standard or test method	Level of immunity test required for professional environment (healthcare facilities)	Compliance level
Electrical fast transient/burst	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV in air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV in air
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz

Table 4 - Test specifications for enclosure port immunity to RF wireless communication devices according to IEC 61000-4-3

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity testing (V/m)
385	380 - 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 - 470	GMRS 460, FRS 460	Frequency modulation ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 - 787	LTE Band 13, 17	217 Hz pulse modulation	0.2	0.3	9
745						
780						
810						
870						
930	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
1720	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	217 Hz pulse modulation	2	0.3	28
1845						
1970						
2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	217 Hz pulse modulation	2	0.3	28
5240	5100 - 5800	WLAN 802.11 a/n	217 Hz pulse modulation	0.2	0.3	9
5500						
5785						

WARNING - Mobile RF communication devices (including antenna cables and external antennas) must not be used closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer, otherwise device performance may be degraded.

4.9 Essential performance

No, the use of the device does not depend on functions necessary to prevent unacceptable risks. The absence or malfunction of electrical stimulation does not pose an unacceptable risk.

5 SERIOUS INCIDENT REPORTING

Any serious incidents occurring in connection with the device must be reported to the manufacturer:

Vygon Italia S.r.l. and to the competent authority of the Member State where the user and/or patient is located.

The list can be accessed on the EU website at the following link:

https://health.ec.europa.eu/medical-devices-sector/new-regulations/contacts_en

"Serious incident" is defined as any incident that, directly or indirectly, has caused, may have caused, or is likely to cause any of the following:

- the death of a patient, user or other person;
- the serious temporary or permanent deterioration in the health of the patient, user or other person;
- a serious threat to public health.

An "incident" is defined as any malfunction or change in the characteristics or performance of a device placed on the market, including misuse due to ergonomic characteristics, as well as any inadequacy in the information provided by the manufacturer and any undesirable side effect.

6 WARRANTY

Vygon Italia S.r.l. guarantees free service (labour and spare parts) for any failure resulting from manufacturing defects for 24 months from the date the device is placed on the market.

Failures caused by drops, impacts, neglect, unsuitable or unauthorised maintenance, or otherwise attributable to improper use of the device are excluded from the warranty.

Parts subject to wear and tear, such as batteries, cables, etc., are also excluded from the warranty.

Maintenance, including non-warranty maintenance, may only be performed by Vygon Italia S.r.l. or authorised centres.

The warranty is excluded in case of interventions carried out by unauthorised third parties. Vygon Italia or its authorised service centre has the right to decide whether to repair or replace certain parts or the entire device. In case of replacement, the replaced parts will be retained by Vygon Italia. A valid supporting document and clear identification of the individual device serial number are required in order to claim warranty repair.

7 DESCRIPTION OF SYMBOLS ON THE LABEL

SYMBOL	TITLE	SYMBOL	TITLE
	CE		Manufacturer
	Medical Device		Catalogue Number
	Date of manufacture		Unique Device ID
	Refer to the instructions for use		Serial Number
	Recycling: electronic equipment		Applied parts, type BF
	Model		

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Vygon Italia S.r.l. e a autoridade competente do Estado Membro no qual o usuário e/ou paciente está localizado.
VYGON ITALIA Srl ▼ **VYGON GROUP**